

TECHSPEC<sup>®</sup> 138mm Telecentric Backlight Illuminator



138mm Telecentric Backlight Illuminator, #35-427

Stock #35-427 3-4 DAYS



1



£1,874.25

ADD TO CART

Qty 1+

£1,874.25

Volume Pricing

[Request Quote](#)

Prices shown are exclusive of VAT/local taxes

**Note:** This item requires accessories for use | [Learn More](#)

Get Product Downloads

**Diameter (mm):**  
160.00

**Length (mm):**  
329.50

**Type of Illumination:**  
Specialty

**Note:**  
Aperture is not adjustable

|                                              |
|----------------------------------------------|
| <b>Manufacturer:</b><br>Edmund Optics®       |
| <b>Geometry:</b><br>Backlight                |
| <b>Beam Diameter (mm):</b><br>138            |
| <b>Mounting Flange Diameter (mm):</b><br>200 |
| <b>Light Mount Inner Diameter (mm):</b><br>8 |

Regulatory Compliance

|                                                            |
|------------------------------------------------------------|
| <b>RoHS:</b><br><a href="#">Compliant</a>                  |
| <b>Certificate of Conformance:</b><br><a href="#">View</a> |

---