

20 x 20mm Glass Checkerboard Target



Glass Checkerboard Calibration Target

Stock #15-974 3-4 DAYS



1



£446.25

ADD TO CART

Qty 1-4

£446.25

Qty 5+

£424.57

Volume Pricing

[Request Quote](#)

⚠ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Coating:

Low Reflective Chrome (reflectivity <25% at 550nm)

Substrate:

Float Glass

Pattern Size (mm):

0.2 x 0.2 ±0.002

Dimensions (mm): 20 x20 ±0.10
Thickness (mm): 1.52 ±0.10
NIST Certification: No

Regulatory Compliance

RoHS 2015: Compliant
Certificate of Conformance: View
Reach 235: Compliant
