

TECHSPEC® 242mm Telecentric Backlight Illuminator



TECHSPEC® Telecentric Backlight Illuminators

Stock #12-229 3-4 DAYS

1 **£4,377.50**

ADD TO CART

Qty 1+

£4,377.50

Volume Pricing

[Request Quote](#)

! Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Diameter (mm):
271.00

Length (mm):
678.00

Type of Illumination:
Specialty

Note:
Aperture is not adjustable

Manufacturer:
Edmund Optics®

Geometry:
Backlight

Beam Diameter (mm):
242

Mounting Flange Diameter (mm):
298.5

Light Mount Inner Diameter (mm):
8

Regulatory Compliance

RoHS:
[Compliant](#)

Certificate of Conformance:
[View](#)
