

TECHSPEC®

400 - 750nm, 25.2 x 35.6mm, Dichroic Longpass Filter Kit (8 Filters)



Stock #15-104

CONTACT US

-

1

+

£1,032.75

ADD TO CART

|           |
|-----------|
| Qty 1+    |
| £1,032.75 |

|                               |
|-------------------------------|
| Volume Pricing                |
| <a href="#">Request Quote</a> |

ⓘ

Prices shown are exclusive of VAT/local taxes

Get Product Downloads



Dimensions (mm):

25.2 x 35.6

Filters Included :

Cut-on Wavelengths: [400nm](#), [450nm](#), [500nm](#), [550nm](#), [600nm](#), [650nm](#), [700nm](#), [750nm](#)

Length (mm):

35.60

Type:

Longpass Filter

|                                                                                                                                                                                  |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Width (mm):</b><br>25.20                                                                                                                                                      |
| <b>Number of Filters:</b><br>8                                                                                                                                                   |
| <b>Environmental Durability:</b><br>Humidity per ML-STD-810H, Section 507.6<br>Temperature per ML-STD-810H, Section 501.7 and 502.7                                              |
| <b>Physical Durability:</b><br>Adhesion per ML-PRF-13830B, Section C.4.5.12<br>Moderate abrasion per ML-PRF-13830B, Section C.4.5.11<br>Cleaning per ML-C-48497A Section 4.5.4.2 |

Regulatory Compliance

|                                                            |
|------------------------------------------------------------|
| <b>RoHS:</b><br><a href="#">Compliant</a>                  |
| <b>Certificate of Conformance:</b><br><a href="#">View</a> |