

TECHSPEC®

800 - 1200nm, 25.2 x 35.6mm, Dichroic Longpass Filter Kit (9 Filters)



Stock #15-105

CONTACT US

-

1

+

£1,126.25

ADD TO CART

Qty 1+
£1,126.25

Volume Pricing
Request Quote

ⓘ

Prices shown are exclusive of VAT/local taxes

Get Product Downloads



Dimensions (mm):

25.2 x 35.6

Filters Included :

Cut-on Wavelengths: [800nm](#), [850nm](#), [900nm](#), [950nm](#), [1000nm](#), [1050nm](#), [1100nm](#), [1150nm](#), [1200nm](#)

Length (mm):

35.60

Type:

Longpass Filter

Width (mm): 25.20
Number of Filters: 9
Environmental Durability: Humidity per ML-STD-810H, Section 507.6 Temperature per ML-STD-810H, Section 501.7 and 502.7
Physical Durability: Adhesion per ML-PRF-13830B, Section C.4.5.12 Moderate abrasion per ML-PRF-13830B, Section C.4.5.11 Cleaning per ML-C-48497A Section 4.5.4.2

Regulatory Compliance

RoHS: Compliant
Certificate of Conformance: View