

TECHSPEC® Mounting Clamp, 2.25" Centerline, 30mm ID



Mbunting Clamp, 2.25" Centerline, #56-871



Stock **#56-871** 3-4 DAYS

- 1 +

£140.25

ADD TO CART

Qty 1+

£140.25

Volume Pricing

[Request Quote](#)

Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Inner Diameter (mm):
30.0

Stock No. of Mounting Clamp:
[#56-870](#) Sold Separately

Type:
Lens Accessory

Regulatory Compliance

RoHS:
[Compliant](#)

Certificate of Conformance:
[View](#)
