

TECHSPEC®

Multi-Element Retaining Ring Set



Multi-Element Retaining Ring Set

Stock **#37-774** 3-4 DAYS

-

1

+

£15.21

ADD TO CART

Qty 1-4

£15.21

Qty 5+

£13.34

Volume Pricing

[Request Quote](#)

⚠ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Function:

Retainer Ring

Thread Type:

M29 x0.75

Construction:

Matte Black Anodized Aluminum (6061-T6)

Thickness (mm): 3.00
Inner Diameter (mm): 25.0

Regulatory Compliance

RoHS 2015: Compliant
Certificate of Conformance: View
Reach 233: Compliant
