


Nikon Tube Lens to C-Mount Holder

See More by [Nikon](#)



Nikon Tube Lens to C-Mount Holder (#11-150)

Stock **#11-150** 3-4 DAYS



1



£126.65

ADD TO CART

Qty 1+

£126.65

Volume Pricing

[Request Quote](#)

Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Length (mm):
37.70

Mounting Threads:
1" x32 TPI (Male and Female)

Type:
Microscope Accessory

Outer Diameter (mm):
50.00

Note:
Requires adapter [#11-148](#) to connect Nikon objective

Regulatory Compliance

RoHS:
[Compliant](#)

Certificate of Conformance:
[View](#)
