

Pedestal Post, English, 50.8mm



Pedestal Post

Stock **#15-848** 3-4 DAYS



1



£22.44

ADD TO CART

Qty 1-4

£22.44

Qty 5+

£20.23

Volume Pricing

[Request Quote](#)

⚠ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Mounting Threads:
8-32

Type:
English

Weight (kg):
0.19

Construction: Stainless Steel
Length (mm): 50.80
Note: With 1/4-20 screw holes on the base
Function: Pedestal Post
Post Diameter (mm): 25.4
Base Diameter (mm): 31.5

Regulatory Compliance

RoHS 2015: Compliant
Certificate of Conformance: View
Reach 233: Compliant