

Pedestal Post, English, 12.7mm



Pedestal Post

Stock #15-845 3-4 DAYS

⊖ 1 ⊕ **£18.70**

ADD TO CART

Qty 1-4

£18.70

Qty 5+

£16.91

Volume Pricing

[Request Quote](#)

ⓘ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Mounting Threads:

8-32

Type:

English

Weight (kg):

0.06

Construction:
Stainless Steel

Length (mm):
12.70

Note:
With 1/4-20 screw holes on the base

Function:
Pedestal Post

Post Diameter (mm):
25.4

Base Diameter (mm):
31.5

Regulatory Compliance

RoHS 2015:
[Compliant](#)

Certificate of Conformance:
[View](#)

Reach 233:
[Compliant](#)
