


TECHSPEC[®] 125mm, Center Drive, Thru-Hole, English Extended Micrometer



TECHSPEC[®] Crossed-Roller Bearing Linear Translation Stages (Thru-Hole Model)



Stock #66-464 5 In Stock

-

1

+

£728⁰⁰

ADD TO CART

Volume Pricing	
Qty 1-2	£728.00 each
Qty 3+	£635.04 each
Need More?	Request Quote

Prices shown are exclusive of VAT/local taxes

Product Downloads

General

Type:

English

Physical & Mechanical Properties

Type of Movement:

Linear (X)

