

1X ZEISS Tube Lens

See More by [ZEISS](#)



1X ZEISS Tube Lens, #13-828

Stock **#13-828** 3-4 DAYS



1



£229.50

ADD TO CART

Qty 1+

£229.50

Volume Pricing

[Request Quote](#)

Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Diameter (mm):
40.00

Focal Length FL (mm):
164.50

Model Number:
425308-0000-000

Length (mm):
12.80

Magnification: 1X
Mounting Threads: M36 x0.75
Type: Tube Lens
Note: Required when using a ZEISS objective with a non-ZEISS microscope
Wavelength Range (nm): 380 - 700
Manufacturer: ZEISS
Field Number (mm): 25
Entrance Pupil Diameter (mm): 31.6

Regulatory Compliance

RoHS: Compliant
Certificate of Conformance: View