

[See all 12 Products in Family](#)

TECHSPEC®

30mm, Metric Rotary Stage



Stock #66-513 5 In Stock

-

1

+

£298⁴⁰

ADD TO CART

Volume Pricing	
Qty 1-2	£298.40 each
Qty 3+	£260.24 each
Need More?	Request Quote

 Prices shown are exclusive of VAT/local taxes

Product Downloads

General

Metric

Type:

Note:

Bottom adapter plate, sold separately, is required to mount 30mm and 40mm stages to a standard breadboard.

Physical & Mechanical Properties

Type of Movement:

