

TECHSPEC® 70mm, 150mm Radius, Metric Goniometer



TECHSPEC® Goniometers



Stock #66-536 **3 In Stock**

⊖ 1 ⊕ £632.⁰⁰

ADD TO CART

Volume Pricing	
Qty 1-2	£632.00 each
Qty 3+	£550.64 each
Need More?	Request Quote

ⓘ Prices shown are exclusive of VAT/local taxes

Product Downloads


General

Type:
Metric

Resolution (arcmin):
5

Physical & Mechanical Properties

70mm Gonimeters



Ø (mm)	English Size	Metric Size	Ø (mm)	Ø (mm)	Ø (mm)
120	M6 S22	M6 S22	130	180	12.5
130	M6 S23	M6 S26	140	130	12.4
160	M6 S24	M6 S27	160	140	12.4