

[See all 4 Products in Family](#)

TECHSPEC®

Aperture Kit for Achromat Prototyping Kit



Aperture Kit for Achromat Prototyping Kit

Stock **#54-720** **1 In Stock**

-

1

+

£134⁴⁰

ADD TO CART

Volume Pricing	
Qty 1-5	£134.40 each
Qty 6-10	£125.04 each
Qty 11-25	£119.28 each
Need More?	Request Quote

i Prices shown are exclusive of VAT/local taxes

Product Downloads

General

Prototyping Kit

Type:

Kit

Function:

Regulatory Compliance

