

TECHSPEC[®] Cell with Two Retainers and Apertures for Achromat Proto Kit



Stock **#54-719** [CONTACT US](#)

-1+

£240^{.00}

ADD TO CART

Volume Pricing	
Qty 1-5	£240.00 each
Qty 6-10	£211.68 each
Qty 11-25	£195.20 each
Need More?	Request Quote

!

 Prices shown are exclusive of VAT/local taxes

Product Downloads

General	
Prototyping Kit	Type:
Kit	Function:

Regulatory Compliance

