

Diffuse reflectance phantom Adipose Tissue 50 x 50 x 10 mm



Stock **#71-280** **2 In Stock**

-

1

+

£940^{.00}

ADD TO CART

Volume Pricing	
Qty 1-4	£940.00 each
Qty 5-9	£846.00 each
Need More?	Request Quote

 Prices shown are exclusive of VAT/local taxes

Product Downloads

General

Type:

Adipose Tissue

Physical & Mechanical Properties

Dimensions (mm):

50 x 50 x 10

Material Properties

Reduced Scattering Coefficient, μ'_s :

7	
0.1	Absorption Coefficient, μ_a :

Regulatory Compliance

View	Certificate of Conformance:
----------------------	-----------------------------

Product Details

- Mimic Muscle, Brain, or Adipose Tissue
- Available in 25 x25 mm and 50 x50 mm Sizes
- Ideal for a Range of Life Science Applications

Diffuse Reflectance Phantoms are solid samples that mimic the optical properties of biological tissues for testing, calibration, and characterization of biophotonics devices. These homogeneous phantoms, with well characterized optical properties at 780nm, simulate the optical properties and mechanical characteristics of biological tissues such as muscle, brain, and adipose tissue. These samples come in two sizes, 25 x 25mm or 50 x 50mm, to accommodate multiple experimental setups. Diffuse Reflectance Phantoms are ideal for spatial frequency domain imaging, diffuse reflectance spectroscopy, diffuse Raman spectroscopy, and biophotonics/biomedical imaging.