

Dalsa Genie Nano 1/4-20 Mounting Adapter

See More by [Teledyne DALSA](#)



Dalsa Genie Nano 1/4-20 Mounting Adapter, #34-966



Stock **#34-966** 3-4 DAYS



1



£42.50

ADD TO CART

Qty 1+

£42.50

Volume Pricing

[Request Quote](#)

Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Type:
Camera Accessory

Regulatory Compliance

RoHS:
[Compliant](#)

REACH 201:
[Compliant](#)

Certificate of Conformance:
[View](#)

