

[See all 8 Products in Family](#)

## 0.5mm Dia., TO-46 Package, InGaAs Photodiode



InGaAs Photodiodes (FC Receptacle , TO-5, TO-46)

Stock **#62-271** [CONTACT US](#)

⊖ 1 ⊕ £84.<sup>00</sup>

**ADD TO CART**

Volume Pricing	
Qty 1-4	£84.80 each
Qty 5-9	£72.00 each
Qty 10-24	£67.60 each
Need More?	<a href="#">Request Quote</a>

ⓘ Prices shown are exclusive of VAT/local taxes

### Product Downloads

#### General

Ball Lens **Type of Optics:**


InGaAs Photodiode **Type:**

**Rise/ Fall Time @  $V_R=5V$  (ns):**

10



Typical InGaAs Detector Response Curve



;