

## Edmund Optics® Lab Storage Kit



Stock #29-092 **4 In Stock**

⊖ 1 ⊕ £72<sup>00</sup>

**ADD TO CART**

Volume Pricing	
Qty 1+	£72.00 each
Need More?	<a href="#">Request Quote</a>

ⓘ Prices shown are exclusive of VAT/local taxes

### Product Downloads

### General

**Contents of Kit:**

- 24 component storage drawers
- 12 optics storage trays
- drawer labels
- marker

### Regulatory Compliance

Certificate of Conformance:

## Product Details

- Maintain an Organized Optics Lab with Wall Mountable Storage
- 24 Component Storage Drawers, 12 Optics Storage Trays, Drawer Labels, and a Marker
- Convenient Storage for Lab Supplies
- Durable Boxes with Labels for Organization

Edmund Optics® Lab Storage Kit provide an ideal solution for protecting and organizing components in an optics lab environment. Each Lab Storage Kit includes 24 component storage drawers, 12 optics storage trays, drawer labels, and a marker. These convenient storage cases can be placed on a lab bench, easily wall mounted, and will help alleviate clutter and prevent loss and damage to precision components

---

;