

[See all 22 Products in Family](#)

## Tube Lens



Stock #18-505 **12 In Stock**

− 1 + £416<sup>00</sup>

**ADD TO CART**

### Volume Pricing

Qty 1-4	£416.00 each
Qty 5+	£395.20 each
Need More?	<a href="#">Request Quote</a>

ⓘ Prices shown are exclusive of VAT/local taxes

### Product Downloads

### General

**Note:**

One required between camera tube and tube lens spacer/filter cube

### Physical & Mechanical Properties

Length (mm):  
16.15

## Optical Properties

120.00 Focal Length FL (mm):

1X Magnification:

## Threading & Mounting

M29 x 0.75 & Ø30mm Dovetail Mounting Threads:

## Regulatory Compliance

[View](#) Certificate of Conformance:

## Product Details


- Up to 43% Smaller than Industry Standard Systems
- Optimized Performance up to 1.1" Sensor
- Ideal for Life Science and Industrial Applications
- Winner of the [2023 Bronze VSD Innovators Award](#)

Edmund Optics® 120i Plan APO Infinity Corrected Objectives are infinite conjugate microscope objectives that offer exceptional image flatness and chromatic correction over 400-700nm. These objectives were designed for a 1.1" sensor to offer a large field of view and reliable image quality for high-throughput imaging applications. Ideal for life science and machine vision applications such as microscopy, flow cytometry, pharmaceutical inspection, and assembly line or fault inspection.

These objectives won the [2023 Bronze VSD Innovators Award](#).

## Technical Information

### Compact Tube Lens System



# Standard Tube Lens System

